

Case study

Analysis of Surfactant - Surface interactions with QSense QCM-D

This case study shows how QSense® QCM-D can be used to study surfactant interactions with surfaces in real-time with nanoscale sensitivity.

Introduction

Surfactants are surface active molecules that modulate the surface tension at air/liquid interfaces as well as phenomena and interactions at liquid/solid interfaces. Surfactants are widely used as important components in detergents and other cleaning agents, in oil processing and recovery, pharmaceutical formations as well as in many other products. It is therefore critical to be able to study the effects of surfactants in real-time and at the nanoscale. The Quartz Crystal Microbalance with Dissipation Monitoring technology (QCM-D) offers this possibility by enabling the analysis of surfactant interactions with surfaces including the kinetics, mass and stability of the surfactant/surface interactions.

Experimental

Two different surfactants were analyzed:

1. Triton X-100: non-ionic detergent with a relatively low Critical Micelle Concentration (CMC) (0.24mM), Mw = 645Da.
2. Octyl-β-D-glucoside (BOG): a non-ionic detergent with a relatively high CMC (20.25mM), Mw = 292Da.

The buffer solution was PBS, pH 7.4. QCM-D sensors with a gold surface were treated with UV ozone prior to the measurement. They were then mounted in the QCM-D flow modules and the modules were filled with buffer solution. Next, the sensors were exposed to a Triton X-100 solution (0.24 mM) and a BOG solution (20 mM), respectively, followed by rinsing with buffer.

Results

The adsorption behavior of surfactants on a gold surface was analyzed in real time. **Figure 1** illustrates the layer thickness as a function of time for the two surfactants. The adsorption of Triton X-100 occurs faster but results in a thinner film than the slower adsorption of BOG which yields a thicker film. The effective mass (thickness) of the formed molecular layers are 280 ng/cm² (25 Å) for Triton X-100 and 740 ng/cm² (68 Å) for BOG, Table 1. Upon rinsing with buffer, a fraction of the bound surfactant de-sorbed.




Figure 1. Effective layer thickness as Triton X-100 and BOG adsorb to the surface. The results show that BOG forms a thicker layer than the Triton X-100. Also the binding rate is different; saturation takes place quicker for the Triton X-100, even if the concentration was much lower

Surfactant	Mass Before Rinse	Mass After Rinse	Film Thickness Before Rinse	Film Thickness After Rinse
Triton X-100	280 ng/cm ²	14 ng/cm ²	2.5 nm	0.1 nm
BOG	740 ng/cm ²	360 ng/cm ²	6.8 nm	3.3 nm

Table 1: Areal mass densities and thicknesses of the adsorbed Triton X-100 and BOG layers before and after rinse.

Particularly for BOG, a considerable amount of surfactant remains on the surface after the rinse

Final masses of 14 ng/cm² for Triton X-100 and 360 ng/cm² for BOG were measured after the rinsing step. This corresponds to average molecular layer thicknesses of 1 Å for Triton X-100 and 33 Å for BOG, Table 1.




Figure 2. The left panel shows the frequency and dissipation responses for the adsorption of Triton X-100, while the right panel shows the responses for the adsorption of BOG. The frequency and energy dissipation responses indicate that the BOG forms a somewhat thicker and much softer layer than the Triton X-100. Overtone number 5, 7 and 9 are displayed (dark to light blue/gold).

Raw Data

In **Figure 2**, the changes in frequency (Δf) and energy dissipation (ΔD) are plotted as a function of time for Triton X (left) and BOG (right). The total shifts in dissipation and spreading of overtones during the adsorption indicates that BOG is resulting in a softer film (likely to have a considerable amount of coupled water) than that obtained with Triton X.

Conclusions

This case study shows that QSense QCM-D technology can detect adsorption of small molecules with $M_w < 1000$ Da. The binding rate, mass and structural properties of the formed molecular layers are analyzed. Conclusions regarding the suitability of

surfactants in different applications can be drawn based on such studies. The surface based QCM-D technology offers the possibility to study interaction processes under a variety of different substrates and experimental conditions.

References

Measurements conducted at Biolin Scientific, Gothenburg, Sweden.

Acknowledgements

Dr Ralf Richter at CIC biomaGUNE, San Sebastian, Spain and Professor Alain Brisson, University of Bordeaux, France are acknowledged for valuable input.

